




SPINEFRONTIER™
WHERE SURGEONS INNOVATE

Arena-L™

Anterior Lumbar Interbody

ALIFT Lumbar IBF System



INDICATIONS:

The SpineFrontier Lumbar Interbody Fusion Device System ALIFT is intended for intervertebral body fusion of the spine of skeletally mature patients, using autogenous bone graft to facilitate fusion. The device is indicated for use in patients with degenerative disc disease (DDD) at one or two contiguous spinal levels from L2-S1. These DDD patients may also have up to Grade I spondylolisthesis or retrolisthesis at the involved level(s).

Degenerative disc disease is defined as a discogenic back pain with degeneration of the disc confirmed by history or radiographic studies. These patients should be skeletally mature and have had six months of non-operative treatment.

The SpineFrontier Lumbar Intervertebral Body Fusion Device System is intended to be used with supplemental spinal fixation system(s) for use in the lumbar spine.

Features and Benefits

- Large bone graft pocket to promote bone fusion
- Can be used with an Anterior or Anteriolateral approach
- Aggressive teeth design prevents expulsion
- Manufactured from PEEK Optima™
- Tapered, beveled design allows for easy insertion
- Low profile, robust inserter



SPINEFRONTIER™
WHERE SURGEONS INNOVATE

500 Cummings Center
Suite 3500
Beverly, MA 01915 U.S.A.

P: + 1 978 232 3990
F: + 1 978 232 3991
innovation@spinefrontier.com

www.spinefrontier.com
DOC81002 rev.01



Learn more at www.les-society.org

Les is enough.