





SPINEFRONTIER™
WHERE SURGEONS INNOVATE



FacetFuse-Allo™

Spine Stabilization System



Features and Benefits

- Surgeon-driven device
- Open or percutaneous minimally invasive technique
- Bone-to-bone fusion through the facet joints.
- No metal in body
- Reduces OR Time

FacetFuse-Allo is a cylindrical shaped implant intended to promote fusion in facet joints in conjunction with an interbody fusion device or supplemental spinal fixation. FacetFuse-Allo is inserted bilaterally through the dorsal side of the lamina, across the facet joint through the pedicle towards the vertebral body.

Processed by:
University of Miami Tissue Bank (UMTB)
1951 NW 7th Avenue, Suite 200, Miami, FL 33136

UMTB is accredited by the American Association of Tissue Banks (AATB)
Distributed by SpineFrontier



SPINEFRONTIER™
WHERE SURGEONS INNOVATE

500 Cummings Center
Suite 3500
Beverly, MA 01915 U.S.A.

P: + 1 978 232 3990
F: + 1 978 232 3991
innovation@spinefrontier.com

www.spinefrontier.com
DOC10083 rev.03



Learn more at www.les-society.org

Les is enough.