




**SPINEFRONTIER™**  
WHERE SURGEONS INNOVATE



# LESR-Tri

Less Exposure Surgical Retractor - Tri-Formation



## Features and Benefits

- Two blades for Cranial/Caudal tissue distraction and a third blade for independent posterior distraction
- Up to 60mm of distraction and radiolucent blades
- Multiple table mount positions
- Robust design for optimum position stability
- Accessories include a bright lighting system, blade stability anchor screws, anatomy matching blade extenders
- Smaller incision than traditional open procedure results in less tissue trauma and scarring
- Can be used with standard fluoroscopic equipment

## Indicated for:

- Procedures requiring posterior, oblique or lateral access to the lumbar and thoracic spine
- LES/MIS procedures



**SPINEFRONTIER™**  
WHERE SURGEONS INNOVATE

500 Cummings Center  
Suite 3500  
Beverly, MA 01915 U.S.A.

P: + 1 978 232 3990  
F: + 1 978 232 3991  
innovation@spinefrontier.com

www.spinefrontier.com  
DOC00705 rev.01



Learn more at [www.les-society.org](http://www.les-society.org)

**Les is enough.**