



SPINEFRONTIER™
WHERE SURGEONS INNOVATE



LESR-X

Less Exposure Surgical Retractor - Cross Formation



Features and Benefits

- Independent blade distraction
- Independent blade angulation up to 45°
- Frame angles to the blades to sit against the patient
- Easy push button release to unlock the blades
- Up to 50mm of distraction
- Radiolucent blades
- Multiple table mount positions
- Robust design for optimum position stability
- Two blades for Cranial/Caudal distraction
- Can be used with standard fluoro equipment or a second retractor for medial/lateral or anterior/posterior retraction
- Accessories include a bright lighting system, blade stability anchor screws, anatomy matching blade extenders
- Smaller incision than traditional open procedure results in less tissue trauma and scarring
- Faster patient recovery with less post-op pain, resulting in shorter hospital stay

Indications

- Procedures requiring posterior, oblique, anterior or lateral access to the lumbar, thoracic or cervical spine
- LES/MIS procedures




SPINEFRONTIER™
WHERE SURGEONS INNOVATE

500 Cummings Center
Suite 3500
Beverly, MA 01915 U.S.A.

P: + 1 978 232 3990
F: + 1 978 232 3991
innovation@spinefrontier.com

www.spinefrontier.com
DOC00704 rev.01



Learn more at www.les-society.org

Les is enough.